


# FLOW CHEMISTRY TEACHING METHOD

## Table of Contents

Chapter 1: Safe Working in a Laboratory	3
1. Safe Working in a Laboratory	4
1.1 Pre-laboratory Preparation	4
1.2 Personal Protective Equipment	4
1.3 General Safety Instructions	4
1.4 Labtrix® Start Safety Instructions	5
1.5 Risk Assessment	5
Chapter 2: Micro Reaction Technology	6
2. Introduction to Micro Reaction Technology	7
2.1 What is a Micro Reactor?	7
2.2 Advantages of Micro Reactors	7
2.3 How are Reactions Performed?	7
2.4 Aim of the Experiments	8
2.5 References	8
Chapter 3: Labtrix® Start, Micro Reaction Screening Apparatus	9
3. General Labtrix® Start Information	10
3.1 Micro Reactors	10
3.2 Reactor Holder and Peripherals	11
3.3 Liquid Handling	13
3.4 Thermal Control	14
3.5 Background Information	14
3.5.1 Reagent Preparation: How to make a Volumetric Stock Solution	14
3.5.2 Loading the Solution into a Gas-tight Syringe	15
3.5.3 Calculation of Reaction Time	15
3.5.4 Qualitative Analysis of Flow Reactions	16
3.5.5 Quantitative Analysis of Flow Reactions	16
3.5.6 Determination of Product Purity	17
Chapter 4: Experimental Section	18
Method 1: Synthesis of 3-Phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol	19
Method 2: Synthesis of ( <i>E</i> )-Ethyl 2-cyano-3-(3,5-dimethoxyphenyl)acrylate	25
Method 3: Synthesis of (3a <i>R</i> ,7a <i>S</i> )-5-Methyl-3a,4,7,7a-tetrahydroisobenzofuran-1,3-dione	31
Method 4: Synthesis of Coumarins under Continuous Flow	37
Method 5: Synthesis of 1-Phenylpropan-1-ol <i>via</i> a Continuous Flow Grignard Reaction	44
APPENDIX I: Labtrix® Start Spare Parts/Part Number List	51
APPENDIX II: Material Compatibilities	52
APPENDIX III: Model Answers	53

## CHAPTER 4: EXPERIMENTAL SECTION

**Method 1. Synthesis of 3-Phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol**

**Introduction:** The aim of the experiment is to optimise the reaction time and temperature to maximise the proportion of 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol, prior to synthesising ~ 100 mg of material for analysis by NMR spectroscopy. Depending on equipment availability, reaction success will be determined using Gas-Chromatography (GC-FID) analysis or Thin Layer Chromatography (TLC).

**1. Risk Assessment**

Students are reminded that they have a legal responsibility to take all necessary precautions to ensure the safety of themselves and others in the Laboratory during the course of this experiment. The substances especially harmful to health for this experiment are;

1. *4,4,4-Trifluoro-1-phenyl-1,3-butanedione*: Harmful by inhalation, in contact with skin and eyes.
2. *Hydrazine monohydrate*: Flammable. May cause cancer. Toxic by inhalation, in contact with skin and if swallowed. May cause sensitization by skin contact. Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment (Exposure limits; STEL = 0.1 ppm).
3. *Phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol*: Harmful by inhalation and in contact with skin.
4. *Ethanol*: Highly flammable liquid and vapour. Harmful by inhalation and ingestion. Delayed effects may include depression of the central nervous system, narcosis and damage to the heart (Exposure limits; TWA = 1000 ppm).
5. *Acetone*: Highly flammable liquid. Harmful by inhalation, ingestion and contact with skin and eyes (Exposure limits; TWA = 500 ppm).
6. *Ethyl acetate*: Highly flammable. Harmful by inhalation; may cause drowsiness (Exposure limits; PEL 400 ppm).
7. *Hexane*: Highly flammable. May cause impaired fertility. Irritating to the eyes, skin and central nervous system.

Please refer to the material safety data sheet for each component to ensure the most up to date safety information is known and adhere to the following precautions;

- Gloves, safety glasses and a laboratory coat must be worn throughout the practical session
- The micro reaction equipment must be set-up within a fume cupboard and used with the sash closed
- The preparation of stock solutions must be conducted within a fume cupboard
- Good laboratory practise should be employed throughout the practical session
- All waste generated from this reaction must be placed in an appropriate labelled waste container

## 2. Reagents and Solvents

Name	Purity (%)	Molecular Weight (g mol <sup>-1</sup> )	Density (g ml <sup>-1</sup> )
4,4,4-Trifluoro-1-phenyl-1,3-butanedione	99	216.16	N/A
Hydrazine monohydrate	+98	50.06	1.027
Ethanol	Laboratory reagent grade	46.06	0.780
Acetone	Laboratory reagent grade	58.08	0.790
Ethyl acetate	Laboratory reagent grade	88.10	0.897
Hexane	Laboratory reagent grade	86.18	0.655
Deuterated dimethyl sulfoxide	99.9	N/A	N/A

## 3. Analysis Conditions

Depending on equipment availability, the micro reactions will be quantified using GC-FID or qualified using TLC under the conditions outlined below.

### 3.1 GC-FID Methodology


Using a gas chromatography (GC-FID) operated under the conditions outlined in [Chapter 3, Section 1.5.5](#), the following retention times are obtained for the solvents, reagents and reaction product;

Name	Retention Time (min)
4,4,4-Trifluoro-1-phenyl-1,3-butanedione	2.82
Ethanol	1.42
3-Phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol	3.52

Upon analysis of the samples generated, calculate the percentage conversion of 4,4,4-trifluoro-1-phenyl-1,3-butanedione to 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol using peak area.

### 3.2 TLC Analysis

Thin layer chromatography is carried out using Kieselgel 60, HF<sub>254</sub> alumina backed TLC plates, with 80:20 hexane:ethyl acetate used as eluent. Visualisation is achieved using fluorescence on exposure to short wave ultra violet light ( $\lambda$  254 nm). The solvent front is marked with a pencil and the R<sub>f</sub> values for the starting material and product. The effect of flow reaction conditions can be assessed qualitatively based on the appearance of the product and consumption of the starting material.



## 4. Experimental Conditions


### 4.1 Stock Solutions

Using the following weights and volumes, prepare the stock solutions volumetrically; see introductory information for details of how to do this-if unsure consult a demonstrator.

Stock No.	Solution	Quantity	Volumetric Flask (ml)	Concentration (M)
1	4,4,4-Trifluoro-1-phenyl-1,3-butanedione	1.08 g	5.00	1.00
2	Hydrazine monohydrate	0.27 ml	5.00	1.10

### 4.2 Micro Reactor Set-up

The reactions will be performed using reactor type 3223 (Volume = 10  $\mu$ l) and the reagents will be introduced into the reactor using the following inlets;



Before starting the reaction;

- Check that the reactor holder contains reactor type 3223; consult a demonstrator if another reactor is installed.
- Fill two 1 ml glass gas-tight syringes with absolute ethanol (EtOH) and attach the luer fittings connected to Inlets 1 and 2. Fill a third 1 ml glass gas-tight syringe with acetone and attach to the luer fitting connected to Inlet 6.
- Mount the two syringes filled with EtOH on Pump 1, the syringe filled with acetone on Pump 2 and set the pump rate of both pumps to  $12.5 \mu\text{l min}^{-1}$ ; collecting the solvent in a waste vial at the reactor outlet for 20 min. Check that solvent is passing through the reactor and back-pressure regulator into the waste vial; if after 2 min no fluid is observed stop both pumps and consult a demonstrator.
- After 20 min, stop the flow on both pumps and remove the syringes. Empty the residual solvent into the appropriate waste container.

Setting up a reaction;

- Fill Syringe 1 with Stock Solution 1 (see above) to the 1 ml mark and attach to Inlet 1 before mounting on Pump 1.
- Fill Syringe 2 with Stock Solution 2 (see above) to the 1 ml mark and attach to Inlet 2 before mounting on Pump 1.
- Fill Syringe 3 with acetone to the 1 ml mark and attach to Inlet 6, mount on Pump 2 and place the outlet tube into a waste vial.
- Set a flow rate of  $12.5 \mu\text{l min}^{-1}$  on Pump 1 and 2 and start both pumps at the same time, collecting the reaction products into a waste vial before proceeding to optimise the reaction conditions (see below).

## 5. Optimising a Continuous Flow Reaction

In order to identify the best reaction conditions to use for the synthesis of 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1*H*-pyrazol-5-ol, the reaction time and temperature must firstly be optimised. To do this, perform the following steps;

- Set the temperature controller to 25 °C and collect the reaction products in a waste vial for 5 min.
- After 5 min, move the outlet tube from the waste vial to a clean labelled vial and collect for 1 min.
- Move the outlet tube to waste and set the flow rate on Pumps 1 and 2 to 5  $\mu\text{l min}^{-1}$  for 5 min.
- After 5 min, move the outlet tube from the waste vial to a clean labelled vial and collect for 2 min.
- Move the outlet tube to waste and set the flow rate on Pumps 1 and 2 to 2.5  $\mu\text{l min}^{-1}$ .
- After 5 min, move the outlet tube from the waste vial to a clean labelled vial and collect for 5 min.
- Move the outlet tube to waste and set the flow rate on Pumps 1 and 2 to 1  $\mu\text{l min}^{-1}$ .
- After 10 min, move the outlet tube from the waste vial to a clean labelled vial and collect for 10 min.

Analyse each sample generated by GC-FID or TLC using Stock Solution 1 as the reference material. Qualitatively determine the effect of reaction time on the formation of 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1*H*-pyrazol-5-ol and record the observation made.

Having evaluated reaction time, perform the following temperature evaluation for a fixed 15 s reaction time;

- Re-fill Syringes 1, 2 and 3 with the respective Stock Solution or acetone and mount on the correct Pump.
- Set Pumps 1 and 2 to 20  $\mu\text{l min}^{-1}$  and collect the reaction products in a waste vial for 2 min.
- Set the temperature controller to 50 °C and collect the reaction products in a waste vial for 5 min.
- After 5 min, move the outlet tube to a clean labelled vial and collect for 1 min, move the outlet tube to waste and set the temperature controller to 100 °C.
- After 5 min, move the outlet tube to a clean labelled vial and collect for 1 min, move the outlet tube to waste and set the temperature controller to 150 °C.
- After 5 min, move the outlet tube to a clean labelled vial and collect for 1 min, move the outlet tube to waste and set the temperature controller to 195 °C.
- After 5 min, move the outlet tube to a clean labelled vial and collect for 1 min, move the outlet tube to waste, leave the temperature controller set to 195 °C and reduce the flow rate on all pumps to 1  $\mu\text{l min}^{-1}$ .

Analyse each sample generated by GC-FID or TLC and qualitatively determine the effect of reactor temperature on the formation of 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1*H*-pyrazol-5-ol and record the observations made.

**DO NOT** touch the micro reactor holder or heater unit until the temperature controller reaches a temperature of 25 °C as there is a risk of injury.

## 6. Determination of Reaction Yield



Under the optimal reaction conditions identified from Section 5, operate the micro reactor for a period of time sufficient to collect 2 ml of reaction mixture; collecting the reaction product into a pre-weighed 10 ml round bottomed flask. Evaporate to dryness on a rotary evaporator, allow the flask to cool and weigh the flask. Determine the mass of 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1*H*-pyrazol-5-ol synthesised. Reduce the set flow rate on each pump to 1  $\mu\text{l min}^{-1}$  and the temperature controller to 25 °C. Once cooled stop the pumps, remove the syringes and empty any residual reagents into the correct waste container. Wash the syringes with acetone, fill to 1 ml and connect to inlets 1, 2 and 6 pumping at 12.5  $\mu\text{l min}^{-1}$  for 20 min to clean the tubing and reactor, in preparation for the next experiment.

## 7. NMR Spectroscopic Analysis

Dissolve 10 mg of the reaction product in  $\text{CDCl}_3$  (+ 2 drops of  $d$ -DMSO) and record a  $^1\text{H}$  and  $^{13}\text{C}$  NMR spectra; where this service is unavailable, interpret the model spectra provided.


### 8. Experimental Write-up

Describe the experimental procedure followed, recording the appearance, weight and percentage yield obtained for the material synthesised and interpret the spectra for the 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol synthesised (if NMR analysis is not available please refer to the model spectra provided). Comment on the differences you would expect in the  $^1\text{H}$  NMR spectra for 4,4,4-trifluoro-1-phenyl-1,3-butanedione in the keto- and enol- forms. Propose a synthetic method for the conversion of 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol to 3-phenyl-5-(trifluoromethyl)-1H-pyrazole.



$^1\text{H}$  NMR (400 MHz) spectra of 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol in  $\text{CDCl}_3$  and  $\text{DMSO-}D_6$ .





<sup>13</sup>C NMR (100 MHz) spectra of 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol in CDCl<sub>3</sub> and DMSO-D<sub>6</sub>.

## APPENDIX III: Model Answers

## Method 1. Synthesis of 3-Phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol

## A. TLC Analysis

Substrate	R <sub>f</sub>
4,4,4-Trifluoro-1-phenyl-1,3-butanedione	0.62
3-Phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol	0.21

## B. NMR Spectroscopy

$\delta_{\text{H}}$  (400 MHz, CDCl<sub>3</sub> and DMSO-D<sub>6</sub>) 3.19 (1H, d, J 17.8, CHH), 3.36 (2H, m, CHH and OH), 7.26 (1H, brs, NH), 7.32-7.42 (2H, m, 2 x ArH), 7.62 (2H, m, 2 x ArH) and 7.88 (1H, s, 1 x ArH).

$\delta_{\text{C}}$  (100 MHz, CDCl<sub>3</sub> and DMSO-D<sub>6</sub>) 41.2 (CH<sub>2</sub>), 90.9 (q, J = 29.8, C<sub>0</sub>OH), 124.2 (q, J = 278.3, CF<sub>3</sub>), 125.6 (2 x CH), 128.6 (3 x CH), 132.3 (C<sub>0</sub>) and 147.4 (C<sub>0</sub>N).

## C. Reaction Yield

Under the optimal conditions of a reaction time of 15 s and a reactor set temperature of 195 °C, 2 ml of collected reaction product will contain 0.153 g of 3-phenyl-5-(trifluoromethyl)-4,5-dihydro-1H-pyrazol-5-ol assuming 100 % conversion (7.673 x 10<sup>-5</sup> g per  $\mu$ l).

## D. Time Taken

The time taken to complete the reaction is estimated to be 5 h.

## Experimental and Equipment Support

For details on system assembly, part numbers and trouble shooting, please refer to the Labtrix Start user manual supplied with the system. For any technical questions relating to the system or the development of additional teaching methods, please contact us at [info@chemtrix.com](mailto:info@chemtrix.com).